Pitch-Patch Tissue Reinforcement

Instructions for use

Hinweise zur Anwendung
Consignes d'utilisation
Istruzioni per l'uso

GMDN 46242

Xiros™ Ltd
Springfield House, Whitehouse Lane,
Leeds, LS19 7UE, UK
Tel: +44 (0) 113 238 7200
Fax: +44 (0) 113 238 7201
enquiries@xiros.co.uk
www.xiros.co.uk
Registered in England No.1664824

All rights reserved. © Xiros™ 2020. Worldwide patents and patents pending. Xiros is a Trademark of Xiros.

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Description

The Pitch-Patch is an implantable device constructed from a knitted fabric with integral eyelets. It is reinforced around its perimeter and around each eyelet to increase the security of suture attachment. It is available in two sizes. It should not be cut to size. It is made from non-absorbable polyester. The device is supplied sterile.

General Use

Use the Pitch-Patch only for the stated Indications listed below and not where Contraindicated.

The Pitch-Patch is available in two sizes for medium and large tears. These are 30 x 20 mm and 35 x 25 mm respectively. It is supplied for single use only.

Refer to the relevant Surgical Technique Manual for any additional Contraindications, Warnings and Precautions that may apply.

Sterility

The Pitch-Patch is supplied sterile and double wrapped. It remains sterile unless the package is damaged or opened. Check the package is undamaged before using the Pitch-Patch. If its sterility is compromised the Pitch-Patch must not be resterilized.

This device has been sterilized using a minimum dose of 25 kGy Gamma Radiation. Do not use after expiration date.

Indications

The Pitch-Patch is a single use device intended to be used for reinforcement of the rotator cuff following or during repair by suture or suture anchors, where weakness exists in the soft tissue.

The Pitch-Patch is not intended to replace normal body structure or provide the full mechanical strength to support the rotator cuff. Sutures, used to repair the tear, and sutures or bone anchors, used to attach the tissue to the bone, provide mechanical strength for the repair. The Pitch-Patch reinforces the rotator cuff and provides a scaffold that is incorporated into the patient's own tissue.

The Pitch-Patch is indicated for patients requiring reinforcement of the rotator cuff where either the tear cannot be completely repaired using normal methods and / or the quality of the soft tissue is poor.

Contraindications

- The Pitch-Patch integrates well with the patient's tissue and consequently may not be suitable where eventual removal of the device is anticipated.
- The Pitch-Patch should not be used if there is any evidence of infection.
- The Pitch-Patch should not be used if there is any evidence of decreased vascularity.
- The Pitch-Patch should not be used if there is any structural or pathological condition of the bone or soft tissue that would be expected to impair healing or secure fixation.
- If the patient is suspected of having any foreign body sensitivity, appropriate tests should be made prior to implantation.
- The use of this device may not be suitable for patients who are skeletally immature. It will not elongate as the patient grows. Bone quality should be assessed and surgery should not be performed on patients with insufficient or immature bone.
- The Pitch-Patch should not be used if the patient is unable or unwilling to restrict activities to prescribed levels or follow the rehabilitation programme during the healing period.
- The Pitch-Patch should not be used in patients with glenohumeral arthritis and inflammatory arthropathy.
- The Pitch-Patch should not be used for any repair that requires rolling, folding, or layering, and which may create a space impermeable to fluid, cells, and blood vessels. Such uses may result in excessive inflammation, drainage, extrusion or infection.

Actions

The Pitch-Patch elicits a minimal acute inflammatory reaction in tissues, followed by a gradual encapsulation of the Pitch-Patch by fibrous connective tissue. The

Pitch-Patch is not absorbed, nor is any significant change in tensile strength known to occur *in vivo*.

Warnings and Precautions

- The general principles of patient selection and sound surgical judgement apply to the reconstruction procedure. Use the device only for the listed Indications and in patients whose conditions are not listed as Contraindications above.
- Care should be taken to select the Pitch-Patch best suited to the size of tear and that it is an appropriate size for both the repair and the available joint space.
- The Pitch-Patch should not be cut to size.
- It is recommended that wherever possible patients are placed on prophylactic antibiotics during anaesthesia, to minimize the risk of latent infections developing at the device site.
- Use sterile technique throughout the procedure.
- Do not resterilize. Discard open, unused Pitch-Patch devices.
- When handling the Pitch-Patch (and appropriate sutures used to repair the tear and attach the Pitch-Patch to the soft tissues), care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to the application of surgical instruments such as forceps or needle holders.
- Users should be familiar with surgical procedures and techniques involving sutures and suture anchors or bone anchors before employing the Pitch-Patch.
- When using transosseous bone tunnels and sutures to help secure the Pitch-Patch, take care to leave an adequate bone thickness between the tunnel and the bone surface or between two parallel tunnels, to provide a sufficient bony bridge to resist expected forces that will be exerted on the bone by the reconstruction. Take into account the quality of the bone.
- Check that any knots tied in the sutures that are used to secure the Pitch-Patch are securely locked.
- Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.
- The patient should be warned not to exceed the prescribed activity levels or to overload the repair before complete healing has occurred.
- The safety and effectiveness of this device has not been established in the following patient populations:
 - Pregnant and breastfeeding women

MRI Safety Information

The Pitch-Patch is MR safe (i.e. an item that poses no known hazards) in all MR environments.

Potential Adverse Effects

- Possible adverse reactions associated with the use of this device include: wound dehiscence, infected wounds, minimal acute inflammatory tissue reaction, and transitory local irritation. If any of these conditions occur, the device may need to be removed at the surgeon's discretion.
- The Pitch-Patch may rupture, or the sutures used to secure the Pitch-Patch may break through the device, or the sutures may break through the soft tissue if prescribed activity levels during rehabilitation are exceeded.

Packaging and Labelling

- The Pitch-Patch implants should be accepted only if the factory packaging and labelling arrive intact.
- Contact the Sales Department if the package has been opened or altered.

Storage

Products must be stored in the original unopened packaging away from moisture, dust, insects, vermin, and extremes of temperature and humidity.

Material Specifications

The Pitch-Patch is made from poly (ethylene terephthalate) fibres (polyester).

Directions for use

- Please refer to the Surgical Technique Manual for full instructions for your chosen procedure. These may be obtained from your local distributor.
- Further information if required may be obtained from your local distributor.

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